

Product datasheet for **MK201474**

Dcp1a Mouse qPCR Template Standard (NM_133761)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_133761
Synonyms:	1110066A22Rik; 4930568L04Rik; AU019772; D14Ertd817e; Mitc1; SMIF
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203381) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)