

Product datasheet for **MK201518**

Folr2 Mouse qPCR Template Standard (NM_008035)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_008035
Synonyms:	FBP; FBP2; Fol; Folb; Folbp-2; Folbp2; FR-; FR-beta; FR-P3
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP204904) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)