

Product datasheet for **MK201532**

Nfkbib Mouse qPCR Template Standard (NM_010908)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010908
Synonyms:	I(Kappa)B(beta); I-kappa-B-beta; Ik; IKapp; IKappaBbeta; IKB; ikB-B; IKB-beta; IkBb; NF-kappa-BIB
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP209063) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)