

Product datasheet for **MK201654**

Vps37a Mouse qPCR Template Standard (NM_033560)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_033560
Synonyms:	2210018P21Rik; 4930592A21Rik; AW261445; D8Erttd531e
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP218158) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)