

Product datasheet for **MK201684**

Mrpl42 Mouse qPCR Template Standard (NM_026065)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026065
Synonyms:	2700009F22Rik; 2900055D03Rik; D10Erttd322e; HSPC204; L31mt; MRP-L31; PTD007; Rpml31
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203766) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)