

Product datasheet for **MK201723**

Bcl10 Mouse qPCR Template Standard (NM_009740)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_009740

Synonyms: A1132454; BCL-10; C81403; CARMEN; cE10; CIPER; CLAP; ME10

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP201253) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)