

Product datasheet for **MK201752**

Clcn5 Mouse qPCR Template Standard (NM_016691)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_016691
Synonyms:	5430408K11Rik; CIC-5; Clc4-1; Clc5; Clcn4-1; D930009B12Rik; DXImx42e; Sfc13; T25545
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP202445) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)