

## Product datasheet for **MK201758**

### Icosl Mouse qPCR Template Standard (NM\_015790)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_015790

**Synonyms:** AU044799; B7-H2; B7h; B7RP-1; BG071784; GI50; GL50; GL50-B; ICOS-L; Icoslg; LICOS; Ly115I

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP206461) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)