

Product datasheet for **MK201810**

Rnf14 Mouse qPCR Template Standard (NM_020012)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_020012
Synonyms:	2310075C09Rik; 2610005D23Rik; AA986456; AU041447; D7Bwg0165e; D18Ert188e; Triad2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP214777) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)