

## Product datasheet for **MK201917**

### Dlgap1 Mouse qPCR Template Standard (NM\_177639)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_177639
Synonyms:	4933422O14Rik; 9630002F18; AI845682; AI848168; BB075781; D17Bwg0511e; Gkap; GKAP/SAPAP
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203716) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)