

## Product datasheet for **MK201981**

### Scn3b Mouse qPCR Template Standard (NM\_153522)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_153522
Synonyms:	1110001K16Rik; 4833414B02Rik; Scnb3
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP219797) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)