

## Product datasheet for **MK202000**

### Igf2bp3 Mouse qPCR Template Standard (NM\_023670)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_023670
Synonyms:	2610101N11Rik; AA522010; AL022933; AU045931; IMP-3; IMP3; Koc13; mimp3; Neilsen
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP206585) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)