

Product datasheet for **MK202010**

Tcf4 Mouse qPCR Template Standard (NM_013685)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_013685
Synonyms:	5730422P05Rik; ASP-I2; bHLHb19; E2-2; E2.2; ITF-2; ITF-2b; ITF2; ME2; MITF-2A; MITF-2B; SEF-2; SEF2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217120) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)