

## Product datasheet for **MK202042**

### Bap1 Mouse qPCR Template Standard (NM\_027088)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_027088
Synonyms:	2300006C11Rik; AA989761; AW553466; mKIAA0272; uch-x4
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201413) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)