

Product datasheet for **MK202044**

Cyb5r3 Mouse qPCR Template Standard (NM_029787)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_029787
Synonyms:	0610016L08Rik; 2500002N19Rik; C85115; Dia-1; Dia1; WU:AL591952.1-001; WU:AL591952.1-002
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203218) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)