

## Product datasheet for **MK202182**

### Spag9 Mouse qPCR Template Standard (NM\_001025430)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001025430
Synonyms:	3110018C07Rik; 4733401I23Rik; 4831406C20Rik; AW552012; Jip4; JLP; JSAP2; JSAP2a
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP216074) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)