

Product datasheet for **MK202258**

Ano6 Mouse qPCR Template Standard (NM_175344)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_175344
Synonyms:	2900059G15Rik; AA407480; AW554778; F730003B03Rik; Tmem16f
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217398) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)