

## Product datasheet for **MK202283**

### Scarb2 Mouse qPCR Template Standard (NM\_007644)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007644
Synonyms:	9330185J12Rik; Cd36l; Cd36l2; LGP; LGP85; LIM; LIMP-2; LIMP II; MLGP85
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP215006) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)