

Product datasheet for **MK202288**

Atp1b1 Mouse qPCR Template Standard (NM_009721)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_009721
Synonyms:	At; Atp4b; Atpb; Atpb-1; NKbeta1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201112) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)