

Product datasheet for **MK202313**

Vta1 Mouse qPCR Template Standard (NM_025418)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_025418

Synonyms: 1110001D18Rik; 1110059P08Rik; AU040813; C85340; LIP5; SBP1

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP218410) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)