

## Product datasheet for **MK202322**

### Coro1c Mouse qPCR Template Standard (NM\_011779)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_011779
Synonyms:	AL022675; AW455561; AW548837
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP202953) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)