

Product datasheet for **MK202347**

Asb13 Mouse qPCR Template Standard (NM_080857)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_080857
Synonyms:	2210015B19Rik; 6430573K02Rik; Asb-13; C85285
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP200708) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)