

Product datasheet for **MK202371**

Rdh11 Mouse qPCR Template Standard (NM_021557)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021557
Synonyms:	2610319N22Rik; AI428145; Arsdr1; AU045252; C85936; CGI-82; HCBP12; M42C60; Mdt1; Psdr1; SCALD
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212723) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)