

Product datasheet for **MK202382**

Mtfp1 Mouse qPCR Template Standard (NM_026443)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026443
Synonyms:	1700020C11Rik; 2610507A21Rik; AV005788; Mtp18
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP213188) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

