

Product datasheet for **MK202408**

H2bc4 Mouse qPCR Template Standard (NM_023422)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_023422
Synonyms:	2610022J01Rik; H2bc6; H2bc8; H2bf; H2bfs; Hist1h2; Hist1h2bc; R74621
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP206195) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)