

Product datasheet for **MK202565**

DbnDD2 Mouse qPCR Template Standard (NM_026797)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_026797

Synonyms: 1110017A21Rik; 2900022J10Rik; AU041050; AW048677; D2Bwg0891e; Nkip; R74724

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP222119) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)