

Product datasheet for **MK202615**

Txn14a Mouse qPCR Template Standard (NM_025299)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_025299

Synonyms: D18Wsu98e; Dim1; ENSMUSG00000057130; Txn14; U5-15kD; U5-15kDa

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP220961) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)