

## Product datasheet for **MK202620**

### Aatf Mouse qPCR Template Standard (NM\_019816)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_019816

**Synonyms:** 4933415H02Rik; 5830465M17Rik; Che-1; Trb

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP200829) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)