

Product datasheet for **MK202641**

B4galnt1 Mouse qPCR Template Standard (NM_008080)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_008080
Synonyms:	4933429D13Rik; Gal-NAc-T; Galgt1; GalNAc-T; GalNAcT; Ggm-2; Ggm2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP219868) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)