

## Product datasheet for **MK202676**

### Adh7 Mouse qPCR Template Standard (NM\_009626)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_009626
Synonyms:	Adh-3; Adh-3e; Adh-3t; Adh3; Adh3-e; Adh3-t; Adh4; Adt-1; AI325182
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP200508) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)