

## Product datasheet for **MK202700**

### **F12 Mouse qPCR Template Standard (NM\_021489)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | NM_021489   |
| <b>Synonyms:</b>        | FX; FXII; HAF   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (MP202638) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| <b>Storage:</b>         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)