

Product datasheet for **MK202725**

Ggt6 Mouse qPCR Template Standard (NM_027819)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_027819
Synonyms:	9030405D14Rik; AI427679; AW490214
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP205232) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)