

Product datasheet for **MK202755**

Atp6v1b1 Mouse qPCR Template Standard (NM_134157)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_134157
Synonyms:	Atp6b1; AW208839; D630003L15; D630030L16Rik; D630039P21Rik; Vpp-3; Vpp3
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201095) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)