

Product datasheet for **MK202814**

PPP3CB Mouse qPCR Template Standard (NM_008914)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_008914

Synonyms: 1110063J16Rik; Calnb; Cnab; CnAbeta

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP212030) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C .



[View online »](#)