

Product datasheet for **MK202871**

Galnt10 Mouse qPCR Template Standard (NM_134189)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_134189
Synonyms:	AU018154; C330012K04Rik; GalNAc-T9; GalNAc-T10; Galnt9; ppGaNTase
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP218037) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)