

Product datasheet for **MK202937**

Sumo1 Mouse qPCR Template Standard (NM_009460)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_009460 |
| Synonyms: | GMP1; PIC1; SENTRIN; SMT3; Smt3C; SMT3H3; SMTP3; SUMO-1; Ubl1 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP215599) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)