

Product datasheet for **MK202941**

Nek8 Mouse qPCR Template Standard (NM_080849)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_080849
Synonyms:	4632401F23Rik; b2b1449C; b2b1449Clo; jck
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP208919) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)