

## Product datasheet for **MK203040**

### **B4galt1 Mouse qPCR Template Standard (NM\_022305)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_022305
<b>Synonyms:</b>	4-GalT; 4-GalT1; B-1; b4Gal-T1; beta-1; beta4Gal-T1; Ga; GalT; Ggtb; Ggtb2
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP218012) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)