

Product datasheet for **MK203060**

Ehd4 Mouse qPCR Template Standard (NM_133838)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_133838
Synonyms:	2210022F10Rik; AI197390; AI846352; AV006278; Past2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP204122) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)