

Product datasheet for **MK203111**

Lrp5 Mouse qPCR Template Standard (NM_008513)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_008513
Synonyms:	BMND1; HBM; LR3; LRP7; mKIAA4142; OPPG
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP207608) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)