

## Product datasheet for **MK203127**

### Hspa8 Mouse qPCR Template Standard (NM\_031165)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | NM_031165   |
| Synonyms:        | 2410008N15Rik; Hsc70; Hsc71; Hsc73; Hsp73; Hspa10   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (MP205962) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)