

Product datasheet for **MK203132**

Epcam Mouse qPCR Template Standard (NM_008532)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_008532
Synonyms:	CD326; EGP; EGP-2; Egp314; Ep-CAM; EpCAM1; GA733-2; gp40; Ly74; Tacsd1; Tacstd1; TROP1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217790) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)