

## Product datasheet for **MK203262**

### **Klf6 Mouse qPCR Template Standard (NM\_011803)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_011803
<b>Synonyms:</b>	AI448727; BCD1; C86813; Copeb; CPBP; FM2; FM6; Ierepo1; Ierepo3; R75280; Zf9
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP207227) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)