

Product datasheet for **MK203280**

Sf3b6 Mouse qPCR Template Standard (NM_025323)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_025323
Synonyms:	0610009D07Rik; 6030419K15Rik; AV001342; Sf3b14
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212962) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)