

Product datasheet for **MK203435**

Dctn4 Mouse qPCR Template Standard (NM_026302)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_026302 |
| Synonyms: | 1110001K06Rik; 4930547K17Rik; C130039E17Rik; p62 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP203981) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)