

Product datasheet for **MK203479**

Cldn10 Mouse qPCR Template Standard (NM_021386)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_021386 |
| Synonyms: | 6720456116Rik; Cldn; Cldn1; Cldn10a; Cldn10b; D14Ertd728; D14Ertd728e |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP220745) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)