

Product datasheet for **MK203669**

Pdlim5 Mouse qPCR Template Standard (NM_019808)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_019808 |
| Synonyms: | 1110001A05Rik; AI987914; C87059; Enh; Enh1; Enh2; Enh3; LIM |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP210649) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)