

Product datasheet for **MK203757**

Mapkap3 Mouse qPCR Template Standard (NM_178907)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_178907

Synonyms: 3PK; AI874665; MAPKAP-K3; MAPKAP3; MapkKapk3; MK-3; MK3

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP208356) mix (1 nmol each primer, sufficient for 200 reactions)

1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA

1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)