

## Product datasheet for **MK203854**

### **Nfatc1 Mouse qPCR Template Standard (NM\_198429)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | NM_198429   |
| <b>Synonyms:</b>        | 2210017P03Rik; AI449492; AV076380; NF-ATc; NFAT2; NFATc; Nfatcb   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (MP220424) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| <b>Storage:</b>         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)