

Product datasheet for **MK203878**

Tap2 Mouse qPCR Template Standard (NM_011530)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_011530
Synonyms:	ABC18; Abcb; Abcb3; AI462429; APT2; Ham; Ham-; Ham-2; Ham2; jas; MTP; MTP2; PS; PSF2; RING11; Tap; Tap-2; Y1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217520) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)